

# Direction for use

## Guardian Leather Care Kit

Guardian Leather Care Kit contains everything that is needed for the maintenance of aniline, semi-aniline and pigmented leather: 250 ml. Guardian Leather Cleaner • 250 ml. Guardian Leather Conditioner • 2 Application Sponges

The products in the Leather Care Kit help to preserve the leather's natural characteristics and appearance. Regular care with Guardian Leather Cleaner and Guardian Leather Conditioner significantly reduces the risk of desiccation, scratching and fading. The products in Guardian Leather Care Kit have been tested by the Danish Technological Institute and have been awarded "The Indoor Climate Label".

### Guardian Leather Conditioner

Anilin, Semianilin and Topcoated leather.

Reduce the risk of scratches and stains by applying Guardian Leather Conditioner to your leather furniture.

Leather conditioner contains substances that will maintain the natural look of the leather.

Use sponge or lint free cloth to massage an even, thin layer with circular movements.

Repeat application 4-6 times a year.

Can be used on most leather surfaces.

Some types of leather can change colour slightly after treatment.

### Guardian Leather Cleaner

Anilin, semi-anilin and Topcoated leather.

Use sponge or lint free cloth to massage leather cleaner into the surface with circular movements.

Apply evenly to avoid discolouration.

It is important to clean the affected area in its entirety.

Cleaning can be repeated after 2 hours.

Cleaning recommended 4-6 times a year.

Treat the leather with Guardian Leather Conditioner after every cleaning and 4-6 times a year.

Can be used on most leather surfaces.

Shake before use.

Keep out of the reach of children.

***Always test the product on a piece that is not visible before use.***

***We hope you enjoy using Guardian Leather Care Kit!***